


DORO®

Cranial Stabilization & Retractor Systems



DORO® – *Simply Better Built ... Simply More Choices.*

„Every product we build represents the commitment of our team to put the needs and benefit of the people we serve first.“



Matthias Schüle
PMI Managing Director

Founded in 1993, pro med instruments (PMI) specialized in developing state-of-the-art cranial stabilization and retractor systems. The DORO® brand soon became a synonym for high quality and unique flexibility worldwide.

With this success comes responsibility. We are committed to continuous improvement of ourselves and our products. That is why we listen to our customers and trade partners. Their visions, wishes and suggestions help to further improve our products, for their goals as well as for the benefit of the patients.

Quality, innovation, continuity, partnership and service are not just words to us. Rather these are lived values that form the foundation of our company's philosophy to provide an extra benefit for surgeons and patients alike.

Today, PMI is present in over 70 countries worldwide and growing. Our privately owned company consists of competent and highly motivated employees, dedicated to improve the quality of our users' and patients' lives.



CONTENTS

Rigid & Non-invasive Headrest Systems

DORO® Aluminum Headrest System	page 4
DORO® Aluminum Adaptors and Base Units	page 6
DORO® Base Unit and Parkbench	page 7
DORO® Easy Connect Navigation Adaptors	page 8
DORO® Teflon® Headrest System	page 9
DORO® Reusable and Disposable Skull Pins	page 10
DORO LUCENT® Disposable Skull Pin	page 11
DORO® Multi-Purpose Skull Clamp	page 12
DORO® Multi-Purpose Systems	page 13
DORO® Horseshoe Headrest	page 14

Intra-Operativ Imaging Solutions

DORO® Radiolucent Horseshoe Headrest	page 14
DORO® Radiolucent Headrest System	page 15
DORO® Intra-Op CT Headrest System	page 17
DORO® Mix & Match	page 19
DORO® Aluminum and Radiolucent Table Adaptors	page 20
DORO® Radiolucent Headrest System Accessories	page 21

Halo & Retractor Systems

DORO LUNA® Retractor System	page 22
DORO® Quick-Clamp J-Arm Retractor System	page 24
DORO® Quick-Clamp Halo Brain Retractor System	page 25
DORO® Quick-Clamp McCue Halo Retractor System	page 26
DORO® Quick-Clamp Retractor System	page 27
DORO® 'Rail-It' and Retractor Accessories	page 28
DORO COBRA™ Flexible Rotary Arms	page 30
DORO® Coated Blades	page 32
DORO® Micro Blades	page 33
DORO® Radiolucent Halo Retractor System	page 34

Accessories

DORO® Accessories	page 35
DORO® Cleaning and Sterilization Trays	page 35

DORO® Aluminum Headrest System



DORO® QR3 Skull Clamp



Skull Pins are not included in the DORO® QR3 Headrest System.
DORO® Skull Pins are shown on page 10-11.

All DORO® Swivel Adaptors
and DORO® Base Units
can be found on page 6.

DORO® QR3 Headrest System (shown with DORO® Skull Pins) Item no. 1001.100 (System)*

Description	Item no.
DORO® QR3 Skull Clamp	1001.001
DORO® Swivel Adaptor Aluminum	3002-00
DORO® Base Unit Aluminum	3001-00

* US System item no. 1001.100-US includes Base Unit Aluminum 3001-00-US.

Key Features

- Highest quality workmanship and material selection.
- Excellent stability and flexibility.
- Quick-Rail® technology: Easily to attach retractor systems and accessories to the skull clamp.
- For supine, prone, lateral and sitting positions.^Δ
- DORO® Mix & Match - interchangeable with DORO® Radiolucent Base Units and Adaptors.^{ΔΔ}
- Navigation Ready for Stryker®, Brainlab and Medtronic.^{ΔΔΔ}

^Δ Sitting positions: Only in conjunction with the DORO® Crossbar Adaptor (shown on page 35).

^Δ Lateral position: In conjunction with the DORO® Parkbench Base Unit (shown on page 7).

^{ΔΔ} Interchangeability is only provided when using the corresponding accessories (shown on page 19).

^{ΔΔΔ} DORO® Easy Connect Navigation Adaptors are shown on page 8.

DORO® QR3 Skull Clamp and 'Rail-It' Accessories



All DORO® Retractor Systems and 'Rail-It' Accessories on page 22-33.



DORO® QR3 Skull Clamp & 'Rail-It' Accessories

Description	Item no.
DORO® QR3 Skull Clamp	1001.001
DORO® QR3 Skull Clamp w/ 2 Quick-Clamps	1001.010
DORO® Quick-Clamp	1202.014
DORO® Quick-Clamp Twin QR Coupling Set	1201.010
DORO® Quick-Clamp IGS Arm Coupling Set	1201.040
DORO® Easy Connect Navigation Adaptor, Stryker®	1204.001
DORO® Easy Connect Navigation Adaptor, Brainlab	1204.002
DORO® Easy Connect Navigation Adaptor, Medtronic	1204.003

Key Features

General:

- 3-pin rigid cranial fixation.
- More stable. The design can bear more load due to state-of-the-art engineering.

Quick-Rail®:

- Clear interface access and larger clamping area for a stronger, solid connection of accessories in the sterile field.
- Easy to locate under multiple layers of sterile drape.
- Saves the surgeon time and effort.
- All 'Rail-It' Accessories are compatible with previous versions of aluminum skull clamps through the Quick-Clamp.

Navigation:

- Built-in starburst, on both sides of the Quick-Rail®, provides direct and secure attachment of movement-sensitive navigation tracking devices.

DORO® Aluminum Adaptors

High quality aluminum adaptors to connect DORO® Aluminum Skull Clamps and Horseshoes to the Base Unit.



DORO® Swivel Adaptor
Item no. 3002-00

- Provides 360° rotation for optimum flexibility in patient positioning.



DORO® Ball Pivot Adaptor
Item no. 3002-50

- Provides 360° rotation and pivoting.
- Provides additional height.
- Unique horizontal inclination up to 90° for maximum flexibility in patient positioning.



DORO® Swivel Adaptor Navigation
Item no. 3002-60

- Provides 360° rotation for optimum flexibility in patient positioning.
- Two additional starbursts for increased options in mounting navigational or similar devices.


DORO® Aluminum Base Units

Largest selection of Base Units to connect head positioning devices to OR tables. Excellent stability and durability.

Easy to install in OR table rod slots.



**DORO® Adjustable Base Unit /
DORO® Adjustable Base Unit US**
Item no. 3001-00/ 3001-00-US



**DORO® Adjustable Base Unit
Muranaka**
Item no. 3001-001



**DORO® Adjustable Base Unit
Mizuho**
Item no. 3001-002



**DORO® Adjustable Base Unit
Eschmann**
Item no. 3001-003



**DORO® Adjustable Base Unit
Australia**
Item no. 3001-008



**DORO® Adjustable Base Unit
Eschmann T-Series**
Item no. 3001-010

DORO® Base Unit Parkbench

Provides additional space for positioning the patient's arm comfortably and safely on the armrest between OR table and Parkbench rod. The Armrest Parkbench is mounted by means of the Universal Side Rail Fitting.



DORO® Base Unit Parkbench
Item no. 3001-006



Armrest Parkbench
Item no. 3001-007



Universal Side Rail Fitting
Item no. 3007-50

Full lateral position may improve hemodynamics and decrease nerve impingement complications.

Key Features

- The patient's arm can be positioned comfortably and safely between OR table and Parkbench.
- For prone, supine and especially useful for lateral positioning.
- Ideal for OR tables without traditional two rod slots.
- Unit may slide under OR table for prone and supine position. Built-in side rails on both arms return the functionality of the original OR table side rails.
- Easy to set-up with one self-locking lever for maximum stability.
- Adjustable in width.

Exclusively from
pmi
pro med instruments



DORO® Easy Connect Navigation Adaptors



DORO® Easy Connect Navigation Adaptor Stryker®

Item no. I204.001

- Connects directly to the QR3 Quick-Rail®.
- For Stryker® Navigation Device.



DORO® Easy Connect Navigation Adaptor Brainlab

Item no. I204.002

- Connects directly to the QR3 Quick-Rail®.
- For Brainlab Navigation Device.



DORO® Easy Connect Navigation Adaptor Medtronic

Item no. I204.003

- Connects directly to the QR3 Quick-Rail®.
- For Medtronic Navigation Device.

DORO® 'Rail-It' Accessories




DORO® Quick-Clamp

Backwards compatible:

The Quick-Clamp can be attached to 1st and 2nd GEN Quick-Rail®.

Allows a secure attachment of the DORO® 'Rail-It' Accessories like the DORO® Easy Connect Navigation Adaptors to the Skull Clamp 3003-00.



DORO® Quick-Rail® Accessories

Starburst and 2nd GEN Quick-Rail® integrated in the Quick-Clamp allows a strong and secure attachment of all DORO® 'Rail-It' and Quick-Rail® Accessories.



IGS Articulated Support Arm

Support arm for instruments supporting image-guided applications.

Connects directly to Quick-Clamp (see page 28).

DORO® QR3 Skull Clamp
- navigation ready.



DORO® Teflon® Headrest System



The Headrest System can be sterilized for high risk infection patients.

DORO® Teflon® Headrest System (shown with DORO® Skull Pins)

Description	Item no.
DORO® Skull Clamp Teflon®	3003-009
DORO® Swivel Adaptor Teflon®	3002-009
DORO® Adjustable Base Unit Teflon®	3001-009



Skull Pins are not included in DORO® Teflon® Headrest System. DORO® Skull Pins are shown on page 10-11.

Key Features

- New QR3 Quick-Rail® technology for quick connection of accessories like the new DORO LUNA® Retractor System as well as Retractor Arms, J-Arms and Halo Retractor Systems.
- Autoclavable.
- Machine washable.
- Conforms to WHO Infection Control Guidelines for cleaning and sterilization.
- Navigation ready.
- 3-pin rigid cranial fixation.
- For supine, prone, lateral or sitting position.^Δ
- Excellent stability and flexibility.
- Exclusive – only available from PMI.

^Δ Sitting positions: Only in conjunction with the DORO® Crossbar Adaptor (shown on page 35).

Teflon is a registered trademark of DuPont de Nemours GmbH.

DORO® Reusable Skull Pins



**DORO® Reusable
Skull Pins
Stainless Steel
Adult**
Color: Black
Item no. 3005-00



**DORO® Reusable
Skull Pins
Stainless Steel
Pediatric**
Color: Black
Item no. 3004-00



**DORO® Reusable
Skull Pins
Titanium
Adult**
Color: Blue
Item no. 3005-50



**DORO® Reusable
Skull Pins
Titanium
Pediatric**
Color: Blue
Item no. 3004-50

Key Features

- Packages of 3 Skull Pins.
- Reusable and durable.
- Precision-manufactured tips.
- Designed for easy installation and removal.

DORO® Disposable Skull Pins



**DORO® Disposable
Skull Pins
Stainless Steel
Adult**
Color: Blue
Item no. 3006-00



**DORO® Disposable
Skull Pins
Stainless Steel
Pediatric**
Color: Yellow
Item no. 3006-10



**DORO® Disposable
Skull Pins
Titanium
Adult**
Color: Turquoise
Item no. 3006-20



**DORO® Disposable
Skull Pins
Titanium
Pediatric**
Color: Green
Item no. 3006-30



**DORO® Disposable
Skull Pins
Stainless Steel
Adult**
Color: Black
Item no. 3006-50



**DORO® Disposable
Skull Pins
Stainless Steel
Pediatric,**
Color: Red
Item no. 3006-60

Key Features

- 12 sterile packs of 3 skull pins (36 pins/box).
- Color-coded for ease of use.
- Convenient packaging with strong protective cap, easy and safe to handle.

Precision-machined skull pin tips for accurate pin placement.

Furrowed grip allows easy removal.

Unique shaft ridges compensates for variation in pin receptacles. Ensures secure fixation of the skull pin.

DORO LUCENT® Disposable Skull Pins



- For intra-operative imaging applications using CT, MRI, X-Ray and Fluoroscopy.
- Designed to reduce imaging artifacts to an absolute minimum.
- The composition of titanium tip and high-performance reinforced polymer guarantees a precise, secure and rigid fixation.
- Sterile packaging.



DORO LUCENT® **Disposable Skull Pins**

Adult, sterile

Color: Black

Item no. 1106.003 (3 pins)

DORO® Multi-Purpose Skull Clamp

Shown with Head Support and Gel Pads



Key Features

- New QR3 Quick-Rail® technology for quick connection of accessories like the new DORO LUNA® Retractor System as well as Retractor Arms, J-Arms and Halo Retractor Systems.
- Provides 3- or 4-pin fixation or non-invasive gel pad cranial positioning.
- Interchangeable with quick-change options: Pin Holders, Gel Ear Rings, Gel Head Rings and adjustable Trio Gel Pad.
- Easy to set-up and use.
- Two torque screw options: 360N/80lbs and 90N/20-lbs.
- Navigation ready.

Configuration A

Standard 3-pin fixation.

Configuration B

4-pin fixation for patients with smaller head sizes where rigid fixation is required.

Configuration C

For all non-invasive applications where cranial stabilization is required. All gel pads are interchangeable.



Also usable with
PMT Halo Adaptors
Item no. 1002.067



Storage Case MPSC System

Item no. 3020-001
(shown with Multi-Purpose
Skull Clamp and Accessories).



Adjustment Key

Item no. 3020-65

DORO®

Multi-Purpose Systems

DORO® Gel Pads



DORO® Gel Pad
Item no. 3020-02



DORO® Trio Gel Pad,
(adjustable)
Item no. 3020-71




DORO® Gel Ear Ring
Open, w/ Base 90 mm
Item no. 3020-68



DORO® Gel Ear Ring
Open, w/ Base 140 mm
Item no. 3020-70



DORO® Gel Head Ring
w/ Base 90 mm
Item no. 3020-61



DORO® Gel Head Ring
w/ Base 140 mm
Item no. 3020-63

DORO® Pin Holder



DORO® Single Pin Holder
Item no. 3020-55



DORO® Dual Pin Holder,
Adult
Item no. 3020-56



DORO® Dual Pin Holder,
Pediatric
Item no. 3020-01



DORO® Dual Pin Holder
Adjustable, Pediatric
Item no. 3020-62

*The Multi-Purpose Skull Clamp Systems are also available with the DORO® Multi-Purpose Skull Clamp, 2nd Generation Quick-Rail®.

DORO® Multi-Purpose Skull Clamp System* Item no. 3020-00 (System)

Description	Item no.
DORO® Disposable Skull Pins, Pediatric (36pcs)	3006-10
DORO® Dual Pin Holder Pediatric	3020-01
DORO® Head Support (w/o Gel Pads)	3020-03
DORO® Multi-Purpose Skull Clamp	3020-50
DORO® Single Pin Holder	3020-55
DORO® Gel Head Ring w/ Base 90 mm (3pcs)	3020-61
DORO® Dual Pin Holder Adjustable, Pediatric	3020-62
DORO® Gel Head Ring w/ Base 140 mm	3020-63
DORO® Adjustment Key	3020-65
DORO® Gel Ear Ring Open w/ Base 90 mm (2pcs)	3020-68
DORO® Trio Gel Pad (adjustable)	3020-71
DORO® Storage Case MPSC System	3020-001

DORO® Abbott Multi-Purpose Skull Clamp System* Pediatric, Item no. 3021-00 (System)

Description	Item no.
DORO® Disposable Skull Pins, Pediatric (36pcs)	3006-10
DORO® Dual Pin Holder Pediatric	3020-01
DORO® Head Support (w/o Gel Pads)	3020-03
DORO® Multi-Purpose Skull Clamp	3020-50
DORO® Extension Bar 90N/20lbs	3020-52
DORO® Single Pin Holder	3020-55
DORO® Gel Head Ring w/ Base 90 mm (3pcs)	3020-61
DORO® Dual Pin Holder Adjustable, Pediatric	3020-62
DORO® Gel Head Ring w/ Base 140 mm	3020-63
DORO® Adjustment Key	3020-65
DORO® Gel Ear Ring Open w/ Base 90 mm (2pcs)	3020-68
DORO® Trio Gel Pad (adjustable)	3020-71
DORO® Storage Case MPSC System	3020-001

DORO® Universal Multi-Purpose Skull Clamp System* Item no. 3022-00 (System)

Description	Item no.
DORO® Disposable Skull Pins, Adult (36pcs)	3006-00
DORO® Disposable Skull Pins, Pediatric (36pcs)	3006-10
DORO® Dual Pin Holder Pediatric	3020-01
DORO® Gel Pad	3020-02
DORO® Head Support (w/o Gel Pads)	3020-03
DORO® Multi-Purpose Skull Clamp	3020-50
DORO® Extension Bar 90N/20lbs	3020-52
DORO® Single Pin Holder	3020-55
DORO® Dual Pin Holder Adult	3020-56
DORO® Gel Head Ring w/ Base 90 mm (3pcs)	3020-61
DORO® Pin Holder Adjustable, Pediatric	3020-62
DORO® Gel Head Ring w/ Base 140 mm	3020-63
DORO® Adjustment Key	3020-65
DORO® Gel Ear Ring Open w/ Base 90 mm (2pcs)	3020-68
DORO® Trio Gel Pad (adjustable)	3020-71
DORO® Storage Case MPSC System	3020-001

DORO® Aluminum Swivel Horseshoe Headrest



DORO® Swivel Horseshoe Headrest, Adult

Item no. 3009-00 w/ Extension Bar
Item no. 3009-02 w/o Extension Bar

- Includes adult size horseshoe gel pads.



DORO® Swivel Horseshoe Headrest, Pediatric

Item no. 3008-00 w/ Extension Bar
Item no. 3008-02 w/o Extension Bar

- Includes pediatric size horseshoe gel pads.

Key Features

- Universal Horseshoe Headrest designed for adult or pediatric applications in prone or supine positioning.
- To be mounted on a DORO® Swivel Adaptor Aluminum or on the DORO® Horseshoe Table Adaptor.
- Gel pads interchangeable for adult or pediatric application.
- Permits vertical and lateral adjustments for proper patient positioning.

Additional Accessories

Description	Item no.
DORO® Horseshoe Pad, Pediatric, right	3008-10
DORO® Horseshoe Pad, Pediatric, left	3008-11
DORO® Horseshoe Pad, Adult, right	3009-10
DORO® Horseshoe Pad, Adult, left	3009-11

DORO® Radiolucent Swivel Horseshoe Headrest

Connects directly to the Swivel Adaptor for 360° rotation for optimum flexibility in patient positioning. Also provides vertical and lateral adjustments.



DORO® Swivel Horseshoe Headrest, Radiolucent, Adult

Item no. 3036-00 w/ Extension Bar

- Includes adult size horseshoe gel pads.



DORO® Swivel Horseshoe Headrest, Radiolucent, Pediatric

Item no. 3035-00 w/ Extension Bar

- Includes pediatric size horseshoe gel pads.

Key Features

- Universal Horseshoe Headrest for adult and pediatric cranial support in the prone or supine position for X-Ray, CT imaging and Fluoroscopy.
- Metal free to minimize artifacts.
- Support carrier features Quick-Rail® for fast connection of retractor arms.
- Interchangeable gel pads for adult and pediatric applications.
- Includes Extension Bar for skeletal traction.

Please see page 18 for more Radiolucent Horseshoe options.

DORO® Radiolucent Headrest System



DORO® Radiolucent Headrest System (shown with DORO® Skull Pins) Item no. 3030-00 (System)*

Description	Item no.
DORO® Disposable Skull Pins Titanium, Adult	3006-20
DORO® Storage Case Headrest System Radiolucent	3030-001
DORO® Base Unit Radiolucent	3031-00
DORO® Transitional Member Radiolucent, long (160 mm)	3032-10
DORO® Transitional Member Radiolucent, short (120 mm)	3032-20
DORO® Swivel Adaptor Radiolucent	3033-00
DORO® Skull Clamp Radiolucent	3034-00

* US System item no. 3030-00-US includes Base Unit Radiolucent 3031-00-US.

Other Accessories

Description	Item no.
DORO® Disposable Skull Pins Titanium, Pediatric (36 pcs)	3006-30
DORO® Transitional Member Radiolucent (to ALU Base Unit)	3032-00
DORO® Swivel Adaptor Radiolucent, short	3033-10
DORO® Swivel Adaptor Radiolucent Interface SM iCT, short	3033-12
Interface Skull Clamp Radiolucent to ALU Swivel Adaptor	3033-50
Interface Skull Clamp Aluminum to RDL Swivel Adaptor	3033-51
Universal DRF Adaptor Radiolucent w/ Key	3038-30

Key Features

- 3-pin cranial fixation for cranial and cervical spine procedures using CT, DSA, X-Ray and Fluoroscopy.
- Skull Clamp, Swivel Adaptors and Transitional Members are completely metal free for improved imaging capabilities.
- Attaches to the side rails of OR tables.
- Compatible and interchangeable with a full line of DORO® Radiolucent and Aluminum Base Units and Accessories (shown on page 17-18).

 For DORO® Skull Pins view page 10-11.

DORO® Radiolucent Skull Clamp



Stryker®
3033-64



Medtronic
3033-65



Brainlab
3033-66



**COMPLETELY METAL FREE
MRI SAFE & COMPATIBLE**

DORO® Skull Clamp Radiolucent & Accessories

Description	Item no.
DORO® Navigation Adaptor, Stryker®	3033-64
DORO® Navigation Adaptor, Medtronic	3033-65
DORO® Navigation Adaptor, Brainlab	3033-66
DORO® Skull Clamp Radiolucent	3034-00
DORO® Flexible Rotary Arm for RDL Quick-Rail® 280 mm	3114-39
DORO® Flexible Rotary Arm for RDL Quick-Rail® 420 mm	3114-49



The DORO® Skull Clamp Radiolucent is shown with DORO® Lucent Skull Pins which are not included. (The DORO® Lucent Skull Pin is shown on page 11).



Key Features

- MRI-SAFE Skull Clamp.
- Highest quality workmanship and material selection.
- Skull Clamp is completely metal free for improved imaging capabilities.
- 3-pin rigid cranial fixation.
- Torque screw has integrated metal free pressure spring.
- Quick-Rail® technology for fast and easy connection of accessories such as Retractor Arms, Retractor Systems and Navigation Adaptors.

DORO® Storage Case Headrest System Radiolucent

Item no. 3030-001 (shown with Headrest System Radiolucent)

DORO® Intra-Op CT Headrest System

Intra-operative imaging enables surgeons to visualize the anatomy of a patient's brain during surgery, compare it to presurgical scans and precisely track the location of their surgical instruments in relation to the anatomy. Particularly useful for neurosurgical procedures such as tumor resections, shunt placements, biopsy procedures, AVM resections, hematoma evacuations, CSF Management and DBS procedures.



DORO® Intra-Op CT Headrest System

Item no. 3030-300 (System)*

Description	Item no.
DORO® Disposable Skull Pins Titanium, Adult	3006-20
DORO® Storage Case Headrest System Radiolucent	3030-001
DORO® Base Unit Radiolucent	3031-00
DORO® Transitional Member Radiolucent, long (160 mm) (2 pcs)	3032-10
DORO® Transitional Member Radiolucent, short (120 mm)	3032-20
DORO® Transitional Member Radiolucent, Bracing	3032-21
DORO® Swivel Adaptor Radiolucent, short	3033-10
DORO® Skull Clamp Radiolucent	3034-00


* US System item no. 3030-300-US includes Base Unit Radiolucent 3031-00-US.

Key Features

- Extra Bracing and Transitional Members allow increased reach for optimum ISO center positioning in scanner.
- Skull Clamp, Swivel Adaptors, Transitional Members and Bracing are completely metal free for improved imaging capabilities.
- Base Unit can be moved below and under the OR table to bring scanner closer to the bed.
- Navigation Tracking Devices can be mounted on the Quick-Rail® of the Skull Clamp or on its starbursts.
- Compatible with most modern OR surgical/imaging tables and with most surgical imaging systems.

 For DORO® Skull Pins view page 10-11.

DORO® Intra-Op CT Horseshoe RS Headrest System



DORO® Intra-Op CT Horseshoe RS Headrest System

Item no. 3036-300 (System)*

Description	Item no.
DORO® Base Unit Radiolucent	3031-00
DORO® Transitional Member Radiolucent, long (160 mm) (2pcs)	3032-10
DORO® Transitional Member Radiolucent, short (120 mm)	3032-20
DORO® Transitional Member Radiolucent, Bracing	3032-21
DORO® Horseshoe RS Headrest RDL, Adult w/ Extension Bar	3036-70

* US System item no. 3036-300-US includes Base Unit Radiolucent 3031-00-US.

Also available for pediatric applications

Description	Item no.
DORO® Horseshoe RS Headrest RDL, Pediatric w/ Extension Bar	3035-70

Accessories for Intra-Op CT and Swivel Horseshoe

Description	Item no.
DORO® Horseshoe Pad, Pediatric, right	3008-10
DORO® Horseshoe Pad, Pediatric, left	3008-11
DORO® Horseshoe Pad, Adult, right	3009-10
DORO® Horseshoe Pad, Adult, left	3009-11

Key Features

- Provides cranial support in the prone or supine position for adults and children for X-Ray, CT imaging and Fluoroscopy.
- Metal free to minimize artifacts.
- Connects directly to the Transitional Member for chin-to-chest positioning.
- Extra Bracing and Transitional Members allow increased reach for optimum ISO center positioning in scanner.
- Base Unit can be moved below and under the OR table to bring patient closer to bed.
- Support carrier features Quick-Rail® for fast connection of retractor arms.
- Permits vertical and lateral adjustments.
- Compatible with pediatric gel pads for younger patients.
- Includes Extension Bar for skeletal traction.



DORO® Mix & Match for Cranial and Spine Applications



CONFIGURATION A

Cervical Spine Applications

Description	Item no.
DORO® QR3 Skull Clamp	1001.001
DORO® Base Unit Radiolucent	3031-00
DORO® Transitional Member Radiolucent, long (160 mm)	3032-10
DORO® Transitional Member Radiolucent, short (120 mm)	3032-20
DORO® Swivel Adaptor Radiolucent	3033-00
DORO® Interface Skull Clamp Aluminum (to RDL Swivel Adaptor)	3033-51



#3033-51 allows the connection of an aluminum Skull Clamp to a radiolucent Swivel Adaptor and Base Unit.



CONFIGURATION B

Cranial and Cervical Spine Applications with CT, DSA, Fluoroscopy

Description	Item no.
DORO® Base Unit Aluminum	3001-00
DORO® Transitional Member Radiolucent, (to ALU Base Unit)	3032-00
DORO® Swivel Adaptor Radiolucent	3033-00
DORO® Skull Clamp Radiolucent	3034-00



#3032-00 allows the connection of a radiolucent Swivel Adaptor and Skull Clamp to an aluminum Base Unit.



CONFIGURATION C

Cranial Applications with CT, DSA, Fluoroscopy

Description	Item no.
DORO® Base Unit Aluminum	3001-00
DORO® Swivel Adaptor Aluminum	3002-00
DORO® Interface Skull Clamp Radiolucent (to ALU Swivel Adaptor)	3033-50
DORO® Skull Clamp Radiolucent	3034-00



#3033-50 allows the connection of a radiolucent Skull Clamp to an aluminum Swivel Adaptor and Base Unit.




DORO® Skull Pins are shown on page 10-11.

DORO® Aluminum OR Table Adaptors




DORO® Table Adaptor

Mounted on square receptacle of OR table.
Item no. 3010-00



DORO® Side Rail OR Table Adaptor

Fits on side rails of OR table lacking inserts.
Item no. 3011-00



DORO® Side Rail OR Table Adaptor

Amsco® Fits 3000, 4000 & 5000 series.
Item no. 3011-10



DORO® Amsco® Headrest Adaptor

Fits 3080 & 3085 series.
Item no. 3011-11

DORO® Radiolucent OR Table Adaptors



DORO® Angiography Table Adaptor **SIEMENS**

for Artis Table
Item no. 4003-00




DORO® Table Adaptor

MAQUET, Radiolucent
for Alphamaquet 1150 and Magnus 1180*
Item no. 3031-50



DORO® Table Adaptor

STERIS, Radiolucent
for STERIS Carbon OR Tables CMAX XRAY
Item no. 4004.002



DORO® Table Adaptor

MEDIFA, Radiolucent
for MEDIFA Hesse OR Table MRT 5600 II
Item no. 4004.001

* Only in conjunction with Maquet adaptor
(Maquet item no. 6005.25A0).

DORO® Radiolucent Headrest System Accessories



**DORO® Swivel Adaptor
Radiolucent**
Item no. 3033-00



**DORO® Swivel Adaptor
Radiolucent, short**
Item no. 3033-10



**DORO® Swivel Adaptor
Radiolucent, Interface SM iCT, short**
Item no. 3033-12



**DORO® Transitional Member
Radiolucent (to ALU Base Unit)**
Item no. 3032-00



**DORO® Transitional Member
Radiolucent, long (160 mm)**
Item no. 3032-10




**DORO® Transitional Member
Radiolucent, short (120 mm)**
Item no. 3032-20



**DORO® Transitional Member
Radiolucent, Bracing**
Item no. 3032-21



**DORO® Interface Skull Clamp
Radiolucent (to ALU Swivel Adaptor)**
Item no. 3033-50



**DORO® Interface Skull Clamp
Aluminum (to RDL Swivel Adaptor)**
Item no. 3033-51

Key



**DORO® Universal DRF Adaptor
Radiolucent (includes Key)**
Item no. 3038-30



**DORO® Safety Guide
Radiolucent**
Item no. 4003.003

** Only in conjunction with DORO® Table Adaptor Maquet, Radiolucent, item no. 3031-50.*

DORO LUNA® Retractor System



DORO
LUNA®



- Quick and easy set-up. Attaches directly to the Quick-Rail®.
- Remarkable flexibility. Can be angled in any desired direction.
- DORO LUNA® Curve can be positioned parallel to the floor independent of skull clamp position.
- Stainless steel yet lightweight.
- Featuring the new DORO COBRA® Flexible Rotary Arms with 360° rotation for precision positioning. Easy and smooth movement. Stable after repositioning under tension.

DORO LUNA® Retractor System

DORO LUNA® Retractor System – Standard System

Item no. I203.100 (System)

Description	Item no.
DORO COBRA® Flexible Rotary Arm 280mm (2 pcs)	I201.006
DORO® Retractor Cleaning and Sterilization Tray	I201.091
DORO LUNA® Quick-Clamp for Quick-Rail® (2 pcs)	I203.001
DORO LUNA® Coupling (4 pcs)	I203.002
DORO LUNA® Rod (2 pcs)	I203.003
DORO LUNA® Curve (2 pcs)	I203.004
DORO LUNA® Cleaning and Sterilization Tray	I203.091
DORO® Blade Set (2 sets)	3116-101

DORO LUNA® Standard Retractor System
& DORO® QR3 Skull Clamp

I203.110

DORO LUNA® Standard Retractor System
& DORO® QR3 Headrest System

I203.120



The DORO LUNA® Retractor System is also available in a bundle with the DORO® QR3 Skull Clamp or the DORO® QR3 Headrest System.

DORO LUNA® Retractor System – Compact System

Item no. I203.200 (System)

Description	Item no.
DORO COBRA® Flexible Rotary Arm 280mm (2 pcs)	I201.006
DORO® Retractor Cleaning and Sterilization Tray	I201.091
DORO LUNA® Quick-Clamp for Quick-Rail®	I203.001
DORO LUNA® Coupling (2 pcs)	I203.002
DORO LUNA® Rod	I203.003
DORO LUNA® Curve	I203.004
DORO LUNA® Cleaning and Sterilization Tray	I203.091
DORO® Blade Set (2 sets)	3116-101

DORO LUNA® Compact Retractor System
& DORO® QR3 Skull Clamp

I203.210

DORO LUNA® Compact Retractor System
& DORO® QR3 Headrest System

I203.220



For more information about the new DORO COBRA® Flexible Rotary Arm please see page 30.




DORO LUNA®
**Quick-Clamp
for Quick-Rail®**
Item no. I203.001



DORO LUNA®
Coupling
Item no. I203.002



DORO LUNA®
Rod
Item no. I203.003



DORO LUNA®
Curve
Item no. I203.004

DORO® Quick-Clamp J-Arm Retractor System

Mounted on DORO® Skull Clamp



Halo Universal Key

Item no. 3110-80



DORO® Quick-Clamp J-Arm Retractor System

Item no. 1203.100T (System)

Description	Item no.
DORO® Quick-Clamp J-Arm Coupling Set (2 pcs)	1203.010
DORO® Halo Universal Key (2 pcs)	3110-80
DORO® Flexible Rotary Arm 280 mm (2 pcs)	3114-37
DORO® Blade Set (5pcs), black coated, 102 mm (2 sets)	3116-101
DORO® Support Handrest J-Arm (2 pcs)	3210-52
DORO® J-Arm 300 mm (2 pcs)	3211-00
DORO® J-Arm Sterilization Case w/ Rack	3214-00

DORO® J-Arm Retractor System w/ C-Clamp

Item no. 3215-00 (System)

Description	Item no.
DORO® Halo Universal Key (2 pcs)	3110-80
DORO® Flexible Arm 230 mm (2 pcs)*	3114-34
DORO® Blade Set (5pcs), black coated, 102 mm (2 sets)	3116-101
DORO® Support Handrest J-Arm (2 pcs)	3210-52
DORO® J-Arm 300 mm (2 pcs)	3211-00
DORO® J-Arm Support Bracket w/ C-Clamp Adaptor (2 pcs)	3213-00
DORO® J-Arm Sterilization Case w/ Rack	3214-00

*Flexible Rotary Arms (3114-37) are also available for this System.



C-Clamp Set-up (3215-00)


Key Features

- One of the most versatile brain retractor systems.
- Attaches directly to the Quick-Rail® of the DORO® Skull Clamp via Quick-Clamps.
- Available in two configurations: Quick-Clamp and C-Clamp version.
- Quick and easy set-up.
- Lightweight yet provides superior stability.
- 1203.100T features flexible rotary arms with full 360° rotation.

Skull Clamp and pins are not included in systems.
More accessories on page 28-33.

DORO® Quick-Clamp Halo Brain Retractor System

Mounted on DORO® Skull Clamp



DORO® Quick-Clamp Halo Retractor System

Item no. 1202.100T (System)


Description	Item no.
DORO® Quick-Clamp Halo Coupling Set (2 pcs)	1202.010
DORO® Halo Ring	3110-01
DORO® Halo Support Rod 8 mm (2 pcs)	3110-08
DORO® Halo Universal Key (2 pcs)	3110-80
DORO® Halo Pattie Tray	3110-85
DORO® Flexible Rotary Arm 280 mm (3 pcs)	3114-37
DORO® Blade Set (5 pcs), black coated, 102 mm (2 sets)	3116-101
DORO® Micro Blade Set (6 pcs)	3116-20
DORO® Halo Sterilization Case w/ Rack	3117-00

DORO® Halo Retractor System w/ C-Clamp

Item no. 3100-00 (System)

Description	Item no.
DORO® Halo Ring	3110-01
DORO® Halo Support Rod 8 mm (2 pcs)	3110-08
DORO® Halo Support Bracket w/ C-Clamp Adaptor (2 pcs)	3110-51
DORO® Halo Universal Key (2 pcs)	3110-80
DORO® Halo Pattie Tray	3110-85
DORO® Flexible Arm 230 mm (3 pcs)*	3114-34
DORO® Blade Set (5 pcs), black coated, 102 mm (2 sets)	3116-101
DORO® Micro Blade Set (6 pcs)	3116-20
DORO® Halo Sterilization Case w/ Rack	3117-00

*Flexible Rotary Arms (3114-37) are also available for this System.



Key Features

- Halo Ring with continuous Quick-Rail® for uninterrupted retractor arm placement.
- Split Rings allow angling for different level of retraction. Half Ring can be removed for greater access.
- Attaches directly to the Quick-Rail® of the DORO® Skull Clamp via Quick-Clamps.
- Available in two configurations: Quick-Clamp and C-Clamp version.
- 1202.100T features flexible rotary arms with full 360° rotation.

Skull Clamp and pins are not included in systems.
More accessories on page 28-33.

DORO® Quick-Clamp McCue Halo Retractor System

Mounted on DORO® Skull Clamp

Halo Universal Key


Item no. 3110-80




Horizontal movement



Horizontal movement



Vertical movement



DORO® Quick-Clamp McCue Halo Retractor System

Item no. 1202.400T (System)

Description	Item no.
DORO® Quick-Clamp Halo Coupling Set (2 pcs)	1202.010
DORO® Half Ring, fixed	3110-02
DORO® Halo Support Rod 8 mm (2 pcs)	3110-08
DORO® Halo Universal Key (2 pcs)	3110-80
DORO® Halo Pattie Tray	3110-85
DORO® Flexible Rotary Arm 280 mm (3 pcs)	3114-37
DORO® Blade Set (5 pcs), black coated, 102 mm (2 sets)	3116-101
DORO® Micro Blade Set (6 pcs)	3116-20
DORO® Halo Sterilization Case w/ Rack	3117-00
DORO® Mini J-Arm, 125 mm (2 pcs)	3209-00
DORO® J-Arm to Halo Adaptor (2 pcs)	3213-001

DORO® McCue J-Arm Halo Retractor System w/ C-Clamp

Item no. 3400-01 (System)

Description	Item no.
DORO® Half Ring, fixed	3110-02
DORO® Halo Support Rod 8 mm (2 pcs)	3110-08
DORO® Halo Support Bracket w/ C-Clamp Adaptor (2 pcs)	3110-51
DORO® Halo Universal Key (2 pcs)	3110-80
DORO® Halo Pattie Tray	3110-85
DORO® Flexible Arm 230 mm (3 pcs)*	3114-34
DORO® Blade Set (5 pcs), black coated, 102 mm (2 sets)	3116-101
DORO® Micro Blade Set (6 pcs)	3116-20
DORO® Halo Sterilization Case w/ Rack	3117-00
DORO® Mini J-Arm, 125 mm (2 pcs)	3209-00
DORO® J-Arm to Halo Adaptor (2 pcs)	3213-001

*Flexible Rotary Arms (3114-37) are also available for this System.


Key Features

- *Unobstructed access to the surgical field.*
An open ring is crucial for skull base applications. Recommended for posterior fossa.
- *Mini J-Arm allows movement vertically and horizontally.*
- *Provides excellent positioning of flexible arms and retractor blades and a better ergonomic work environment for surgeons.*
- *Quick-Rail® on the Halo and J-Arm provides uninterrupted access for retractor arm placement.*
- *Attaches directly to the Quick-Rail® of the DORO® Skull Clamp via Quick-Clamps.*
- *Available in two configurations: Quick-Clamp and C-Clamp version.*
- *1202.400T features flexible rotary arms with full 360° rotation.*

Skull Clamp and pins are not included in systems.
More accessories on page 28-33.

DORO® Quick-Clamp Retractor System

Mounted on DORO® Skull Clamp



Halo Universal Key

Item no. 3110-80



DORO® Quick-Clamp Retractor System

Item no. 1201.100T (System)

Description	Item no.
DORO® Quick-Clamp Twin QR Coupling Set	1201.010
DORO® Quick-Clamp	1202.014
DORO® Halo Universal Key (2pcs)	3110-80
DORO® Halo Pattie Tray	3110-85
DORO® Flexible Rotary Arm 350 mm (2pcs)	3114-38
DORO® Flexible Rotary Arm 420 mm	3114-48
DORO® Quick-Rail® Storage and Sterilization Case	3114-58
DORO® Blade Set (5pcs), black coated, 102 mm (2sets)	3116-101

Other Key Accessories


Description	Item no.
DORO® QR3 Skull Clamp	1001.001
DORO® Multi-Purpose Skull Clamp	3020-50
DORO® Flexible Rotary Arm 200 mm	3114-36
DORO® Flexible Rotary Arm 280 mm	3114-37
DORO® Blade Set (5pcs) conical, black coated, 152 mm	3116-10

Key Features

- Attaches directly to the Quick-Rail® of the DORO® Skull Clamp via Quick-Clamps.
- For quick, secure connection and repositioning of retractor arms.
- Quick-Clamp provides superior stability even over several layers of sterile drape.
- 1201.100T features flexible rotary arms with full 360° rotation.

Skull Clamp and pins are not included in systems.
More accessories on page 28-33.

DORO® 'Rail-It' and Retractor Accessories



Key Features

- Larger clamping area provides a strong, solid connection.
- Maintains superior stability of accessories even over multiple layers of sterile drape.
- Built-in starburst and 12mm Quick-Rail®.
- Made of stainless steel. Can be sterilized.



**DORO® Quick-Clamp
Twin QR Coupling Set**

Item no. I201.010



**DORO® Quick-Clamp
Coupling Twin QR**

Item no. I201.043



**DORO® Quick-Clamp
Halo Coupling Set**

Item no. I202.010



DORO® Quick-Clamp

Item no. I202.014



**DORO® Quick-Clamp
Coupling Universal**

Item no. I202.015



**DORO® Halo Support
Rod Holder**

Item no. I202.016



**DORO® Quick-Clamp J-Arm
Coupling Set**

Item no. I203.010



**DORO® Quick-Clamp
Coupling J-Arm**

Item no. I203.015




**DORO® Quick-Clamp
IGS Arm Coupling Set**

Item no. I201.040



**DORO® Quick-Clamp IGS
Articulated Support Arm**

Item no. I201.041



**DORO® Quick-Clamp w/
Interface Mizuho-Type Retractors**

Item no. I201.045



**Halo Support Bracket w/
C-Clamp Adaptor**

Item no. 3110-51

DORO® 'Rail-It' and Retractor Accessories



DORO® Halo Ring
Item no. 3110-01



DORO® Half Ring, fixed
Item no. 3110-02



DORO® Halo Support Rod 8mm
Item no. 3110-08



DORO® J-Arm 300 mm
Item no. 3211-00



DORO® Mini J-Arm 125 mm
Item no. 3209-00



DORO® Support Handrest J-Arm
Item no. 3210-52



**J-Arm Support Bracket w/
C-Clamp Adaptor**
Item no. 3213-00



DORO® J-Arm to Halo Adaptor
Item no. 3213-001



Clamp Body for Skin Hooks
Item no. 3114-53



Halo Pattie Tray
Item no. 3110-85



Mini Vice
Item no. 3114-54



Halo Universal Key
Item no. 3110-80

- Smooth and very precise movement.
- Stable after repositioning under tension.
- „Clicks“ on DORO LUNA® Curve. Easy to use and move to the desired location.
- Provides 360° body rotation while mounted for extra flexibility in positioning.
- Thumb lock - no tool necessary.
- Minimized obstruction in the surgical field.


DORO COBRA® Flexible Rotary Arms

Description	Item no.
DORO COBRA® Flexible Rotary Arm 200 mm	I201.005
DORO COBRA® Flexible Rotary Arm 280 mm	I201.006
DORO COBRA® Flexible Rotary Arm 350 mm	I201.007
DORO COBRA® Flexible Rotary Arm 420 mm	I201.008




DORO
COBRA²®

DORO® Flexible Arms




DORO® Flexible Arms

Description	Item no.
DORO® Flexible Arm 150 mm	3114-33
DORO® Flexible Arm 230 mm	3114-34
DORO® Flexible Arm 300 mm	3114-35




DORO® Flexible Rotary Arms

Description	Item no.
DORO® Flexible Rotary Arm 200 mm	3114-36
DORO® Flexible Rotary Arm 280 mm	3114-37
DORO® Flexible Rotary Arm 350 mm	3114-38
DORO® Flexible Rotary Arm 420 mm (not shown)	3114-48



DORO® Flexible Rotary Arms for RDL Quick-Rail®

Description	Item no.
DORO® Flexible Rotary Arm for RDL Quick-Rail® 280 mm	3114-39
DORO® Flexible Rotary Arm for RDL Quick-Rail® 420 mm	3114-49




DORO® Flexible Arms w/ interface for other brands' halos


Description	Item no.
DORO® Flexible Arm 150 mm	3114-03
DORO® Flexible Arm 230 mm	3114-04
DORO® Flexible Arm 300 mm	3114-05

Key Features

- Standard and flexible rotary arms for optimal positioning and unlimited retraction within the surgical field.
- Easy and fast to connect and reposition on the Quick-Rail®.
- Maintain superior stability even over multiple layers of sterile drape through Quick-Clamp.
- Flexible rotary arms provide full 360° rotation.
- Compatible with all DORO® McCue Halo, Halo, J-Arm and Brain Retractor Systems.
- Connects directly to all DORO® Skull Clamps via Quick-Clamp.



DORO® Coated Blades




DORO® Blade Set, black coated

Description	Item no.
Set 5pcs/ 152 mm	3116-10
152 x 6 mm	3116-11
152 x 10 mm	3116-12
152 x 16 mm	3116-13
152 x 19 mm	3116-14
152 x 25 mm	3116-15
152 x 13 mm (not in set)	3116-16

DORO® Blade Set, black coated

Description	Item no.
Set 5pcs/ 102 mm	3116-101
102 x 6 mm	3116-111
102 x 10 mm	3116-121
102 x 16 mm	3116-131
102 x 19 mm	3116-141
102 x 25 mm	3116-151
102 x 13 mm (not in set)	3116-161




DORO® Blade Set, conical, black coated

Description	Item no.
Set 5pcs/ 152 mm	3116-60
152 x 2 mm	3116-61
152 x 3 mm	3116-62
152 x 4 mm	3116-63
152 x 5 mm	3116-64
152 x 6 mm	3116-65

DORO® Blade Set, conical, black coated

Description	Item no.
Set 5pcs/ 102 mm	3116-601
102 x 2 mm	3116-611
102 x 3 mm	3116-621
102 x 4 mm	3116-631
102 x 5 mm	3116-641
102 x 6 mm	3116-651




DORO® Blade Set, rounded tip, black coated

Description	Item no.
Set 5pcs/ 102 mm	3116-70
102 x 6 mm	3116-71
102 x 10 mm	3116-72
102 x 16 mm	3116-73
102 x 19 mm	3116-74
102 x 25 mm	3116-75

Key Features

- Extensive accessories program with black mat finish retractor blades, in conical, rounded and square configurations.





DORO® Micro Blade Set

Description	Item no.
Set 6pcs	3116-20
114 x 2 mm	3116-21
114 x 4 mm	3116-22
140 x 2 mm	3116-23
140 x 4 mm	3116-24
165 x 2 mm	3116-25
165 x 4 mm	3116-26




DORO® Micro Blade Set, serrated

Description	Item no.
Set 6pcs	3116-30
114 x 2 mm	3116-31
114 x 4 mm	3116-32
140 x 2 mm	3116-33
140 x 4 mm	3116-34
165 x 2 mm	3116-35
165 x 4 mm	3116-36




DORO® Micro Blade Set, hooked

Description	Item no.
Set 6pcs	3116-40
114 x 2 mm	3116-41
114 x 4 mm	3116-42
140 x 2 mm	3116-43
140 x 4 mm	3116-44
165 x 2 mm	3116-45
165 x 4 mm	3116-46



DORO® Micro Blade Set, tapered, length 127 mm


Description	Item no.
Set 4pcs	3116-50
3,2 x 9,5 mm	3116-51
5,1 x 15,9 mm	3116-52
6,3 x 19,0 mm	3116-53
7,9 x 25,4 mm	3116-54



DORO® Adaptor for Flat Brain Spatula

Item no. 1201.046

Also fits retractor arms shown on page 31.



DORO® Radiolucent Halo Retractor System

Mounted on DORO® Skull Clamp Radiolucent

DORO® Halo Retractor System Radiolucent Item no. 3037.00 (System)

Description	Item no.
DORO® Halo Ring Radiolucent	3037-01
DORO® Rod and Fixation Assembly, Radiolucent (right & left rod)	3037-04
DORO® Halo Universal Key (2 pcs)	3110-80
DORO® Halo Pattie Tray	3110-85
DORO® Flexible Rotary Arm for RDL Quick-Rail® 280 mm (3 pcs)	3114-39
Blade Set (5 pcs), black coated (2 sets)	3116-101



Halo Quick-Rail® Adaptor
RDL/Halo Alu (3114-59)

Quick-Rail® Adaptor J-Arm RDL/Alu
(3114-60) (not shown)



DORO® Open Halo Retractor System Radiolucent Item no. 3038-00 (System)

Description	Item no.
DORO® Quarter Ring Radiolucent, variable	3037-02
DORO® Half Ring Radiolucent, fixed	3037-03
DORO® Rod and Fixation Assembly, Radiolucent (right & left rod)	3037-04
DORO® Halo Universal Key (2 pcs)	3110-80
DORO® Halo Pattie Tray	3110-85
DORO® Flexible Rotary Arm for RDL Quick-Rail® 280 mm (3 pcs)	3114-39
DORO® Blade Set (5 pcs), black coated (2 sets)	3116-101



DORO® Radiolucent
Open Halo
Retractor System



DORO® Halo Ring
Radiolucent
Item no. 3037-01



DORO® Quarter Ring
Radiolucent
Item no. 3037-02
DORO® Half Ring
Radiolucent
Item no. 3037-03

Key Features

- Mounts directly to the DORO® radiolucent Skull Clamp via Quick-Rail®.
- Halo is completely metal free and may remain in place during intra-operative imaging.
- Split Ring allows angling for different level of retraction. Half ring can be removed for greater access.
- The DORO® Open Halo Retractor System features a half-ring with a quarter ring for maximum unobstructed access to the surgical field. Quarter ring can be removed for greater access.

DORO® Accessories



DORO® Crossbar Adaptor

Item no. 3007-00

For sitting positions.



Universal Side Rail Fitting

Item no. 3007-50

For attachment of 3007-00 or 3001-007.



Disposable Headrest Protective Drape

Item no. 3010-10

Non-sterile. Fits all Headrest Systems.



Headrest Wall Storage Unit

Item no. 3013-00

32" w x 37" h.



Cervical Spine Support

Item no. 3012-00

Neck support for anterior spine surgery.



DORO® Table Adaptor Horseshoe

Item no. 3009-16


- Connects the DORO® Horseshoe Headrest directly to the OR table without Swivel Adaptor and Base Unit. Adjustable in multiple planes.
- Also available with a hexagon interface. DORO® Table Adaptor Horseshoe hexagon. Item no. 3009-17

For more accessories visit our website
www.pmisurgical.com

DORO® Cleaning and Sterilization Trays

DORO® Cleaning and Sterilization Trays

Description	Item no.
DORO® Retractor Cleaning and Sterilization Tray	1201.091
DORO LUNA® Cleaning and Sterilization Tray	1203.091
Quick-Rail® Storage and Sterilization Case	3114-58
Halo Sterilization Case w/ Rack (shown)	3117-00
J-Arm Sterilization Case w/ Rack	3214-00
Halo Sterilization Case (RDL System)	6300-13



www.pmisurgical.com



DORO® Representative

Headquarters Germany

pro med instruments GmbH | Bötzingen Straße 38 | 79111 Freiburg, Germany
Phone +49 761 384 222 10 | E-mail pmi@pmisurgical.com

Customer Service & Warehouse Germany

pro med instruments GmbH | Waldmatten 6 | 79224 Umkirch, Germany
Phone +49 761 384 222 11 | E-mail pmi@pmisurgical.com

US Subsidiary

pro med instruments, Inc. | 4529 S.E. 16th Place, Suite 101 | Cape Coral, FL 33904
Toll Free 877 225 4086 | E-mail pmi.us@pmisurgical.com